



Sending your Itero scans to OrthoDenco Labs

1. Open a new Rx and enter patient's information:

Treatment Information

Case

* Patient First Name:

* Patient Last Name:

Chart #:

* Case Type: ▼

Dentist



* Name:

* License #:

Signature:

Scanning Order


* Scanning Order: ▼

6043 Kimberly Blvd, Suite N • North Luaderdale, Fl 33068 • 1-800-315-8829

www.orthodenco.com

2. Select **iRecord** under Case Type:



Treatment Information

Case

* Patient First Name: John

* Patient Last Name: Smith

Chart #:

* Case Type: [Dropdown Menu]

Scanning Order

* Scanning Order: [Dropdown Menu]

Dentist

* Name: Doctor, Orthodontist


* License #: 40987

Signature:

Buttons: Scan Now, Cancel, Scan Later

Note: After selecting iRecord, the ship to section will become visible.

3. Select the lab from the drop-down:



Treatment Information

Case

* Patient First Name: John

* Patient Last Name: Smith

Chart #:

* Case Type: iRecord

Scanning Order

* Scanning Order: Lower Jaw First

Dentist

* Name: Doctor, Orthodontist

* License #: 40987

Signature:

Ship To

Lab: Orthodenco labs

To send to a lab, select the lab in "Ship To"

Scan Now Cancel Scan Later

After scanning the patient, send case as usual.

Please keep in mind:

- If no lab is selected, you will be receiving a standard iRecord file.
- If the **incorrect** lab is selected, the case cannot be redirected.
- If no lab has been related to the office, the **"Ship To"** field will be grayed out and cannot be used.
- If OrthoDenco is not in your ship to list, please call Itero to 800-577-8767 and ask rep to add OrthoDenco to your profile, the lab's ID is 14193.

6043 Kimberly Blvd, Suite N • North Luaderdale, Fl 33068 • 1-800-315-8829

www.orthodenco.com