Procedure Prescription Form

OrthoDenco Orthodontic Lab FLDL # 10742 www.orthodenco.com

Toll Free: 1.800.315.8829 - Direct: 954.917.0499 - Fax: 954.935.9547


Practice Info:

Patient:

Date:

DUE BY (Before appointment)

Case Notes:



Signature and Lic#

Ship Date:
WB ST MT AP
Duter
FI PO S&R



Sustainable Appliance Solutions

Materials Used in Restoration:

- -Retainers: Acrylic monomer & polymer, stainless steel wire and/or screw, silver solder
- -All metal: Stainless steel wire, screw, bands and/or crowns, silver solder
- -Functionals: Stainless steel wire and/or acrylic monomer and polymer
- -Hard Night Guard: Acrylic monomer and polymer and stainless steel wire
- -Soft Night Guard: Ethylen Venylacetat, Vinyl